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YS0119 immFast™ Incomplete Freund's Adjuvant

YS0120 immFast™ Complete Freund's Adjuvant

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Adjuvants are nonspecific stimulators used to improved immune response in conjunction with an antigen. Adjuvants are mixed and injected with an antigen to prevent catabolism and help increase the immune response by localizing the antigen for an extended time and attracting the appropriate cells (T cells, B cells and APC) to interact with it.

Freund's adjuvants produce a stronger, longer lasting immunogenic response compared to other adjuvants. They are easy-to-use water-in-oil emulsions and are supplied in two forms. **immFast™ Freund's Complete Adjuvant** is the form that contains killed cells of *Mycobacterium butyricum* to enhance the immune response. The form that does not contain this bacterium is known as **immFast™ Incomplete Freund's Adjuvant**. The complete adjuvant is used in initial injections, and the incomplete form is used to boost.

【Storage】

Upon receipt store at 2-8°C. Product is shipped at ambient temperature.

【Procedure】

1. Mix equal volumes of immunogen and adjuvant.
2. Mix by forcing the adjuvant-immunogen mixture through a small orifice. Push the mixture back and forth between the syringes for approximately 10-15 minutes until a thick emulsion develops. The mixing is complete if a drop of the emulsion does not disperse on the surface of water. If it does disperse, continue mixing until this test is successful. When this procedure is complete, the sample is ready for injection.

【Attentions 注意事项】

- Products are supplied as immunogen-grade emulsions.
- Avoid direct contact with this adjuvant.
- These reagents are for use in laboratory animals only.
- Injection at multiple sites is generally favorable.
- To avoid anaphylaxis, adjuvants are not generally used for intravenous injection.